

## Anti-Disks large homolog 1 Antibody (6H872)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	6H872
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of SAP97 expression in SH-SY5Y cell lysate.
Application:	WB
Recommended	WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Disks large homolog 1
Synonyms:	SAP97;DLG1;hDlg;DLGH1

### Research Background

Essential multidomain scaffolding protein required for normal development. Recruits channels, receptors and signaling molecules to discrete plasma membrane domains in polarized cells.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481