

Anti-DNAI1 Antibody (8B226)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	8B226
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of Dynein intermediate chain 1 expression in mouse testis lysate.
Application:	WB
Recommended	WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human DNAI1
Synonyms:	ICS;ICS1;Immotile cilia syndrome 1;CILD1;PCD;DNAI 1

Research Background

Part of the dynein complex of respiratory cilia.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481