

Anti-DRP 3 Antibody (1F568)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	1F568
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of CRMP4 expression in Human fetal brain lysate.
Application:	ICC/IF,IP,WB
Recommended	ICC/IF: 1:50-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human DRP 3
Synonyms:	ULIP1;LCRMP;DPYSL3;CRMP4;DRP3;TUC4

Research Background

Necessary for signaling by class 3 semaphorins and subsequent remodeling of the cytoskeleton. Plays a role in axon guidance, neuronal growth cone collapse and cell migration.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481