

Anti-DRP 5 Antibody (1X44)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	1X44
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of CRMP5 expression in SH-SY5Y cell lysate.
Application:	IHC,IP,WB
Recommended	IHC: 1:100-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human DRP 5
Synonyms:	CRAM;Ulip6;DRP5;DPYSL5

Research Background

Involved in the negative regulation of dendrite outgrowth.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481