

Anti-Dynactin 5 Antibody (2E328)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	2E328
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of DCTN5 expression in MCF7 cell lysate.
Application:	ICC/IF,IHC,IHC-Fr,IP,WB
Recommended	ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Dynactin 5
Synonyms:	Dynactin 5 (p25);DCTN 5;Dynactin 4;DCTN5

Research Background

Part of the dynactin complex that activates the molecular motor dynein for ultra-processive transport along microtubules.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481