

Anti-eEF1B gamma Antibody (9B195)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	9B195
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of eEF1G expression in (1) 293T cell lysate; (2) HepG2 cell lysate; (3) NIH/3T3 cell lysate.
Application:	FCM,ICC/IF,IHC,IHC-Fr,IP,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human eEF1B gamma
Synonyms:	EF 1G;eEF1By;eEF 1B gamma;EF1 gamma;EEF1G;EF1G;GIG35;eEF1B γ

Research Background

Probably plays a role in anchoring the complex to other cellular components.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481