

Anti-eIF2B5 Antibody (3Y807)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3Y807
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of eIF2B epsilon expression in Hela cell lysate.
Application:	ICC/IF,IHC,IHC-Fr,WB
Recommended	ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human eIF2B5
Synonyms:	CACH;CLE;EIF2BE;EIF 2B;EIF2Bepsilon;EIF2B;EIF2B5

Research Background

Catalyzes the exchange of eukaryotic initiation factor 2-bound GDP for GTP.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481