

Anti-Enkurin Antibody (8J398)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Clone:	8J398
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of C10orf63 expression in Mouse testis lysate.
Application:	WB
Recommended	WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Enkurin
Synonyms:	C10orf63;ENKUR

Research Background

Adapter that functions to localize a calcium-sensitive signal transduction machinery in sperm to a calcium-permeable ion channel. Microtubule inner protein (MIP) part of the dynein-decorated doublet microtubules (DMTs) in cilia axoneme, which is required for motile cilia beating.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481