

Anti-ERp27 Antibody (2P688)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	2P688
Purification:	Affinity-chromatography

Applications

Application:	ICC/IF,IHC,IHC-Fr,WB
Recommended	ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human ERp27
Synonyms:	C12orf46;PDIA8;ERp 27

Research Background

Specifically binds unfolded proteins and may recruit protein disulfide isomerase PDIA3 to unfolded substrates. Binds protein substrates via a hydrophobic pocket in the C-terminal domain.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481