

Anti-EXOG Antibody (6K257)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6K257
Purification:	Affinity-chromatography

Applications

Application:	IP,WB
Recommended	IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human EXOG
Synonyms:	ENDOGL2;Endo G like 1;ENGL;Nuclease EXOG;ENGLB;ENGLA;Endonuclease G like 1; Endonuclease G like 2;Exog;ENDOGL1

Research Background

Endo/exonuclease with nicking activity towards supercoiled DNA, a preference for single-stranded DNA and 5'-3' exonuclease activity.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481