

Anti-FKBP8 Antibody (3I904)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3I904
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of FKBP38 expression in Jurkat cell lysate.
Application:	FCM,ICC/IF,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human FKBP8
Synonyms:	hFKBP38;FKBP 38;Sam11;Rotamase;FKBPR38;Fkbp8

Research Background

Constitutively inactive PPIase, which becomes active when bound to calmodulin and calcium. Seems to act as a chaperone for BCL2, targets it to the mitochondria and modulates its phosphorylation state.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481