

Anti-FYN binding protein 1 Antibody (5Z436)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Clone:	5Z436
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of FYB expression in TF-1 cell lysate.
Application:	ICC/IF,IHC,IHC-Fr,WB
Recommended	ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human FYN binding protein 1
Synonyms:	ADAP;Fyb;SLAP130;p120/p130;FYN binding protein

Research Background

Acts as an adapter protein of the FYN and LCP2 signaling cascades in T-cells. May play a role in linking T-cell signaling to remodeling of the actin cytoskeleton.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481