

Anti-GCS alpha 1 Antibody (8Z864)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	8Z864
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of GCYA3 expression in Human fetal kidney lysate.
Application:	WB
Recommended	WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human GCS alpha 1
Synonyms:	GUCSA3;GUC1A3;GCS α 1;GCS α 1;GC-S-alpha-1;GCSa;GCSa1;GCS α 1;GCS α 1;GCS-alpha-3;GUCY1A3;GCS a1;GC-SA3;GUCA3;GCS a 1;GUCY1A1

Research Background

Belongs to the adenylyl cyclase class-4/guanylyl cyclase family. Contains 1 guanylate cyclase domain.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481