

## Anti-Glycodelin Antibody (4Q810)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4Q810
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of Placental Protein 14 / PAEP expression in TPA-treated K562 cell lysate.
Application:	WB
Recommended	WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Glycodelin
Synonyms:	Glycodelin A;PAEG;PEG;PAEP;GdA;Placental protein 14;GdF;Glycodelin S;Glycodelin F;PP14;GdS;PEP

### Research Background

Glycoprotein that regulates critical steps during fertilization and also has immunomodulatory effects. Four glycoforms, namely glycodelin-S, -A, -F and -C have been identified in reproductive tissues that differ in glycosylation and biological activity. Glycodelin-A has contraceptive and immunosuppressive activities.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481