

Anti-Glycoprotein Ibalpha Antibody (9G377)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	9G377
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of CD42b expression in human spleen lysate.
Application:	IHC,WB
Recommended	IHC: 1:100-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Glycoprotein Ibalpha
Synonyms:	VWDP;CD42b alpha;GP Ib alpha;Antigen CD42b alpha;DBPLT3;BDPLT3;GP1bA;GP1BA;GP1b alpha;BDPLT1;GP1B;BSS

Research Background

GP-Ib, a surface membrane protein of platelets, participates in the formation of platelet plugs by binding to the A1 domain of vWF, which is already bound to the subendothelium.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481