

## Anti-Golgin subfamily A member 1 Antibody (1C527)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	1C527
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of Golgin 97 expression in HeLa cell lysate.
Application:	IP,WB
Recommended	IP: 1:20-1:50; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Golgin subfamily A member 1
Synonyms:	Gap junction protein alpha 4 37kD; GOLGA1; Golgin97; Golgin subfamily A member 1; Golgi autoantigen golgin subfamily a 1; Golgin A1

### Research Background

Involved in vesicular trafficking at the Golgi apparatus level. Involved in endosome-to-Golgi trafficking.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481