

Anti-HAUS7 Antibody (9W107)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	9W107
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of UCHL5IP expression in Jurkat cell lysate.
Application:	IHC,IP,WB
Recommended	IHC: 1:100-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human HAUS7
Synonyms:	UCHL5IP;UIP1;HSXQ28ORF;Haus7;UCH37 interacting protein 1;HAUS 7

Research Background

Contributes to mitotic spindle assembly, maintenance of centrosome integrity and completion of cytokinesis as part of the HAUS augmin-like complex.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481