

Anti-HNF 1 beta Antibody (2Q284)

Product Details

| | |
|---------------|-------------------------|
| Ig Type: | IgG |
| Reactivity: | Human |
| Conjugation: | Unconjugated |
| Clone: | 2Q284 |
| Purification: | Affinity-chromatography |

Applications

| | |
|--------------------|---|
| Verified Activity: | Western blot analysis of HNF1 beta expression in SK-OV-3 cell lysate. |
| Application: | WB |
| Recommended | WB: 1:1000-1:2000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------|---|
| Immunogen: | A synthesized peptide: human HNF 1 beta |
| Synonyms: | HNF2;TCF 2;HNF 1β;HNF1 beta;HNF 1 β;HNF1β;TCF2;LFB3;MODY5;HNF1 b;HNF 1 b;vHNF1;HPC11;HNF1 β;HNF 2;HNF1B;FJHN;HNF 1B |

Research Background

Transcription factor, probably binds to the inverted palindrome 5'-GTTAATNATTAAC-3'.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481