

Anti-Hsf2 Antibody (4F201)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4F201
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of HSF2 expression in HeLa cell lysate.
Application:	FCM,ICC/IF,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Hsf2
Synonyms:	HSTF 2;HSF2;Heat shock factor protein 2;Hsf 2

Research Background

DNA-binding protein that specifically binds heat shock promoter elements (HSE) and activates transcription. In higher eukaryotes, HSF is unable to bind to the HSE unless the cells are heat shocked.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481