

Anti-Hsp70 12A Antibody (5E590)

Product Details

| | |
|---------------|-------------------------|
| Ig Type: | IgG |
| Reactivity: | Human,Mouse,Rat |
| Conjugation: | Unconjugated |
| Clone: | 5E590 |
| Purification: | Affinity-chromatography |

Applications

| | |
|--------------------|--|
| Verified Activity: | Western blot analysis of HSPA12A expression in U-87 MG cell lysate. |
| Application: | ICC/IF,IHC,IHC-Fr,IP,WB |
| Recommended | ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------|--|
| Immunogen: | A synthesized peptide: human Hsp70 12A |
| Synonyms: | FLJ13874;heat shock 70kD protein 12A;Heat shock 70 kDa protein 12A;heat shock 70kDa protein 12A;KIAA0417;Hspa12a;HS12A |

Research Background

Adapter protein for SORL1, but not SORT1. Delays SORL1 internalization and affects SORL1 subcellular localization.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481