

## Anti-Human Ig lambda light chain Antibody (7Q744)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7Q744
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of Lambda Light chain expression in human plasma lysate.
Application:	FCM,ICC/IF,IHC,IHC-Fr,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Human Ig lambda light chain
Synonyms:	IGLC2;BJP;IGLC1;Constant region of lambda light chains;Paraprotein;IGLC6;Immunoglobulin lambda constant 1;IGLV;IGLC3;Mcg marker;Bence Jones protein;Ig lambda chain C regions;IGLC

### Research Background

Constant region of immunoglobulin light chains. Immunoglobulins, also known as antibodies, are membrane-bound or secreted glycoproteins produced by B lymphocytes.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481