

Anti-IL-5 receptor alpha Antibody (3Z313)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3Z313
Purification:	Affinity-chromatography

Applications

Application:	WB
Recommended	WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human IL 5 receptor alpha
Synonyms:	CD125;CDw125;HSIL5R3;IL5R;IL5RA

Research Background

Cell surface receptor that plays an important role in the survival, differentiation, and chemotaxis of eosinophils.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481