

Anti-Jagged 2 Antibody (8A590)

Product Details

| | |
|---------------|-------------------------|
| Ig Type: | IgG |
| Reactivity: | Human,Mouse |
| Conjugation: | Unconjugated |
| Clone: | 8A590 |
| Purification: | Affinity-chromatography |

Applications

| | |
|--------------------|--|
| Verified Activity: | Western blot analysis of JAG2 expression in HUVEC cell lysate. |
| Application: | WB |
| Recommended | WB: 1:1000-1:2000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------|---------------------------------------|
| Immunogen: | A synthesized peptide: human Jagged 2 |
| Synonyms: | Jagged2;Syndactylism;SER2;HJ2;JAG2 |

Research Background

Putative Notch ligand involved in the mediation of Notch signaling. Involved in limb development.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481