

## Anti-Kinetochore protein NDC80 homolog Antibody (8G505)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	8G505
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of HEC1 expression in Jurkat cell lysate.
Application:	IP,WB
Recommended	IP: 1:20-1:50; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Kinetochore protein NDC80 homolog
Synonyms:	TID3;HsHec1;hsNDC80;ndc80;Kinetochore protein Hec1;HEC;KNTC2;Kinetochore-associated protein 2

### Research Background

Acts as a component of the essential kinetochore-associated NDC80 complex, which is required for chromosome segregation and spindle checkpoint activity.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481