

Anti-Laminin gamma 1 Antibody (60814)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	60814
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of Laminin gamma 1 expression in A431 cell lysate.
Application:	FCM,WB
Recommended	FCM: 1:20-1:100; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Laminin gamma 1 LAMB2;Laminin-6 subunit gamma;lanc1;Laminin B2 chain;Laminin subunit gamma-1;
Synonyms:	Laminin-3 subunit gamma;Laminin-1 subunit gamma;Laminin-11 subunit gamma;Laminin-9 subunit gamma;Laminin-8 subunit gamma;Laminin-2 subunit gamma;Laminin-7 subunit gamma;Laminin-4 subunit gamma;Laminin-10 subunit gamma

Research Background

Binding to cells via a high affinity receptor, laminin is thought to mediate the attachment, migration and organization of cells into tissues during embryonic development by interacting with other extracellular matrix components.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481