

## Anti-Microsomal GST 1 Antibody (90923)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	90923
Purification:	Affinity-chromatography

### Applications

Verified Activity:	All lanes use the Microsomal GST 1 Antibody, mAb, Rabbit at 1:2000 dilution for 1 hour at room temperature.
Application:	FCM,ICC/IF,IHC,IHC-Fr,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Microsomal GST 1
Synonyms:	MGST 1;GST12;MGST I;MGST;MGST1

### Research Background

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481