

## Anti-MRRF Antibody (1D20)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Clone:	1D20
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of MRRF expression in HeLa cell lysate.
Application:	IHC,WB
Recommended	IHC: 1:100-1:200; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human MRRF
Synonyms:	mitochondrial;RRF;MTRRF;Ribosome recycling factor

### Research Background

Responsible for the disassembly of ribosomes from messenger RNA at the termination of mitochondrial protein biosynthesis. Acts in collaboration with GFM2. Promotes mitochondrial ribosome recycling by dissolution of intersubunit contacts.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481