

Anti-Mucin 5B Antibody (1A277)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	1A277
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of MUC5B expression in Hela cell lysate.
Application:	FCM,WB
Recommended	FCM: 1:20-1:100; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Mucin 5B
Synonyms:	MUC5;MG1;MUC5B;MUC9;Mucin5B

Research Background

Gel-forming mucin that is thought to contribute to the lubricating and viscoelastic properties of whole saliva and cervical mucus.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481