

Anti-Myosin 19 Antibody (20132)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	20132
Purification:	Affinity-chromatography

Applications

Application:	ICC/IF,WB
Recommended	ICC/IF: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Myosin 19
Synonyms:	myosin XIX;Myohd1;myo19

Research Background

Actin-based motor molecule with ATPase activity that localizes to the mitochondrion outer membrane. Motor protein that moves towards the plus-end of actin filaments. Required for mitochondrial inheritance during mitosis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481