

Anti-Myosin 1B Antibody (7X205)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	7X205
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of MYO1B expression in (1) 293 cell lysate; (2) Mouse heart lysate.
Application:	IHC,WB
Recommended	IHC: 1:100-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Myosin 1B
Synonyms:	MYH 1c;Myo1b;MM1a;myr1;MMI alpha

Research Background

Motor protein that may participate in process critical to neuronal development and function such as cell migration, neurite outgrowth and vesicular transport.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481