

## Anti-NDK B Antibody (7Q435)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	7Q435
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of NME2 expression in LnCaP cell lysate.
Application:	FCM,ICC/IF,IHC,IHC-Fr,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human NDK B
Synonyms:	C myc purine binding transcription factor PUF;NDPKB;NM23B;PUF;Nucleoside diphosphate kinase B;MGC2212;NDKB;NDP kinase B

### Research Background

NM23 nucleoside diphosphate kinase A. Major role in the synthesis of nucleoside triphosphates other than ATP. Holoenzyme is a hexamer of the A and B subunits. Mutations have been identified in aggressive neuroblastomas. Two transcript variants encoding different isoforms have been described.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481