

## Anti-ODF2 Antibody (7U913)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7U913
Purification:	Affinity-chromatography

### Applications

Verified Activity:	All lanes use the ODF2 Antibody, mAb, Rabbit at 1:2000 dilution for 1 hour at room temperature.
Application:	FCM,WB
Recommended	FCM: 1:20-1:100; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human ODF2
Synonyms:	Cenexin 1;ODF 84;ODF84;ODF 2;Cenexin;ODF2/1;odf2;KKT4;ODF2/2;CT134

### Research Background

Seems to be a major component of sperm tail outer dense fibers (ODF). ODFs are filamentous structures located on the outside of the axoneme in the midpiece and principal piece of the mammalian sperm tail and may help to maintain the passive elastic structures and elastic recoil of the sperm tail.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481