

Anti-p90 S6 kinase alpha 6 Antibody (4S739)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	4S739
Purification:	Affinity-chromatography

Applications

Verified Activity:	All lanes use the p90 S6 kinase alpha 6 Antibody, mAb, Rabbit at 1:1000 dilution for 1 hour at room temperature.
Application:	FCM,ICC/IF,IHC,IHC-Fr,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human p90 S6 kinase alpha 6
Synonyms:	pp90RSK4;p90RSK6;rps6ka6;Ribosomal S6 kinase 4;RSK 4

Research Background

Constitutively active serine/threonine-protein kinase that exhibits growth-factor-independent kinase activity and that may participate in p53/TP53-dependent cell growth arrest signaling and play an inhibitory role during embryogenesis.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481