

## Anti-Peroxin 11B Antibody (3F657)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	3F657
Purification:	Affinity-chromatography

### Applications

Application:	IHC,WB
Recommended	IHC: 1:100-1:200; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Peroxin 11B
Synonyms:	PEX11 beta;PEX11 $\beta$ ;PEX11 $\beta$ ;Pex11p beta;PEX14B;Peroxin11B;PEX11B

### Research Background

Involved in peroxisomal proliferation. May regulate peroxisomes division by recruiting the dynamin-related GTPase DNM1L to the peroxisomal membrane.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481