

## Anti-PPIF Antibody (9N512)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	9N512
Purification:	Affinity-chromatography

### Applications

Verified Activity:	All lanes use the PPIF Antibody, mAb, Rabbit at 1:1000 dilution for 1 hour at room temperature.
Application:	IHC,IP,WB
Recommended	IHC: 1:100-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human PPIF
Synonyms:	Rotamase;Ppif;Mitochondrial cyclophilin;cyclophilin D;hCyp3;CYP3;PPIase F;Cyclophilin F;PPIase;Cyp M;Cyclophilin 3;Cyp D;CypD

### Research Background

PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481