

## Anti-PROS1 Antibody (2B892)

### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	2B892
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of PROS1 expression in Human plasma lysate.
Application:	WB
Recommended	WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human PROS1
Synonyms:	PSA;PS24;THPH6;PROS;PS22;PS21;PS23;Preproprotein S;proS1;THPH5;PS25;PS26;PROS 1

### Research Background

Anticoagulant plasma protein; it is a cofactor to activated protein C in the degradation of coagulation factors Va and VIIIa. It helps to prevent coagulation and stimulating fibrinolysis.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481