

Anti-PSMC4 Antibody (9S443)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	9S443
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of Tbp7 expression in HeLa cell lysate.
Application:	ICC/IF,WB
Recommended	ICC/IF: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human PSMC4
Synonyms:	MIP224;PSMC 4;RPT3;TBP 7;Tat binding protein 7

Research Background

The 26S protease is involved in the ATP-dependent degradation of ubiquitinated proteins. The regulatory (or ATPase) complex confers ATP dependency and substrate specificity to the 26S complex.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481