

Anti-RASEF Antibody (5U91)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	5U91
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of RAB45 expression in (1) A431 cell lysate; (2) Raw264.7 cell lysate.
Application:	FCM,ICC/IF,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human RASEF
Synonyms:	RAB45

Research Background

Binds predominantly GDP, and also GTP. Acts as a dynein adapter protein that activates dynein-mediated transport and dynein-dynactin motility on microtubules.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481