

## Anti-Reticulon 4 receptor Antibody (8P746)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Clone:	8P746
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of Nogo Receptor expression in HepG2 cell lysate.
Application:	IP,WB
Recommended	IP: 1:20-1:50; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Reticulon 4 receptor
Synonyms:	Nogo 66 receptor;NgR;Nogo receptor;Rtn4r;Nogor

### Research Background

Receptor for RTN4, OMG and MAG,Functions as a receptor for the sialylated gangliosides GT1b and GM1.Besides, functions as a receptor for chondroitin sulfate proteoglycans.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481