

Anti-SART1 Antibody (1B989)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	1B989
Purification:	Affinity-chromatography

Applications

Application:	FCM,ICC/IF,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; WB: 1:1000-1:5000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human SART1 (hSART-1);SART 1;U4/U6.U5 tri-snRNP-associated 110 kDa protein;Squamous cell carcinoma antigen recognized by T-cells 1 (SART-1;U4/U6.U5 tri-snRNP-associated protein 1;SNU66 homolog (hSnu66)
Synonyms:	

Research Background

Plays a role in mRNA splicing as a component of the U4/U6-U5 tri-snRNP, one of the building blocks of the spliceosome. May also bind to DNA.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481