

## Anti-Secretogranin 1 Antibody (70952)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	70952
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of Chromogranin B expression in (1) HeLa cell lysate; (2) Raw 264.7 cell lysate.
Application:	WB
Recommended	WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Secretogranin 1
Synonyms:	ChromograninB;Secretogranin B;SCG1;CgB;Sgl;Chgb;Secretogranin1;CCB peptide

### Research Background

Secretogranin-1 is a neuroendocrine secretory granule protein, which may be the precursor for other biologically active peptides.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481