

Anti-SLC12A1 Antibody (2X576)

Product Details

| | |
|---------------|-------------------------|
| Ig Type: | IgG |
| Reactivity: | Human |
| Conjugation: | Unconjugated |
| Clone: | 2X576 |
| Purification: | Affinity-chromatography |

Applications

| | |
|--------------|-------------------|
| Application: | WB |
| Recommended | WB: 1:1000-1:2000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------|--------------------------------------|
| Immunogen: | A synthesized peptide: human SLC12A1 |
| Synonyms: | Slc12a1;NKCC2;BSC1 |

Research Background

Renal sodium, potassium and chloride ion cotransporter that mediates the transepithelial NaCl reabsorption in the thick ascending limb and plays an essential role in the urinary concentration and volume regulation.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481