

Anti-SLC4A1 Antibody (5J772)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	5J772
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of Band 3 / CD233 expression in Red blood cell lysate.
Application:	ICC/IF,IHC,IHC-Fr,WB
Recommended	ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human SLC4A1
Synonyms:	FR;SLC4A1;CD233;BND3;DI;RTA1A;SW;WD;EPB3;WD1;AE1;EMPB3;WR

Research Background

Functions both as a transporter that mediates electroneutral anion exchange across the cell membrane and as a structural protein.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481