

## Anti-SLP 2 Antibody (6M66)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	6M66
Purification:	Affinity-chromatography

### Applications

Application:	FCM,ICC/IF,IP,WB
Recommended	FCM: 1:20-1:100; ICC/IF: 1:50-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human SLP 2
Synonyms:	Stomatin like 2;STOML2;Paratarg 7;Stomatin like protein 2;EPB72 like protein 2;SLP2;HSPC108; Stomatin-like protein 2;Stomatin (EPB72) like 2;Paraprotein target 7

### Research Background

Mitochondrial protein that probably regulates the biogenesis and the activity of mitochondria. Stimulates cardiolipin biosynthesis, binds cardiolipin-enriched membranes where it recruits and stabilizes some proteins including prohibitin and may therefore act in the organization of functional microdomains in mitochondrial membranes.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481