

Anti-Sorting nexin 4 Antibody (8W66)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	8W66
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of SNX4 expression in (1) A431 cell lysate; (2) Raw264.7 cell lysate.
Application:	IP,WB
Recommended	IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Sorting nexin 4
Synonyms:	snx4;ATG24B

Research Background

Involved in the regulation of endocytosis and in several stages of intracellular trafficking. Plays a role in recycling endocytosed transferrin receptor and prevent its degradation. Involved in autophagosome assembly by regulating trafficking and recycling of phospholipid scramblase ATG9A.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481