

Anti-SS5R Antibody (8P692)

Product Details

| | |
|---------------|-------------------------|
| Ig Type: | IgG |
| Reactivity: | Human,Mouse,Rat |
| Conjugation: | Unconjugated |
| Clone: | 8P692 |
| Purification: | Affinity-chromatography |

Applications

| | |
|--------------------|---|
| Verified Activity: | All lanes use the SS5R Antibody, mAb, Rabbit at 1:2000 dilution for 1 hour at room temperature. |
| Application: | FCM,WB |
| Recommended | FCM: 1:20-1:100; WB: 1:1000-1:2000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------|-----------------------------------|
| Immunogen: | A synthesized peptide: human SS5R |
| Synonyms: | SSTR5;SS5 R |

Research Background

Receptor for somatostatin 28 and to a lesser extent for somatostatin-14. The activity of this receptor is mediated by G proteins which inhibit adenylyl cyclase.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481