

Anti-Stanniocalcin 1 Antibody (6C47)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Clone:	6C47
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of STC1 expression in Human Prostate lysate.
Application:	IP,WB
Recommended	IP: 1:20-1:50; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Stanniocalcin 1
Synonyms:	Stc1;STC-1;STC;Stanniocalcin

Research Background

Stimulates renal phosphate reabsorption, and could therefore prevent hypercalcemia.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481