

## Anti-Tat SF1 Antibody (6F224)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Clone:	6F224
Purification:	Affinity-chromatography

### Applications

Verified Activity:	All lanes use the Tat SF1 Antibody, mAb, Rabbit at 1:1000 dilution for 1 hour at room temperature.
Application:	ICC/IF,IHC,IHC-Fr,IP,WB
Recommended	ICC/IF: 1:50-1:200; IHC: 1:100-1:200; IHC-Fr: 1:50-1:200; IP: 1:20-1:50; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Tat SF1
Synonyms:	TAT SF1;HTSF1;HTATSF1

### Research Background

Component of the 175 U2 SnRNP complex of the spliceosome, a large ribonucleoprotein complex that removes introns from transcribed pre-mRNAs.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481