

Anti-Tax1 binding protein 1 Antibody (9Q237)

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	9Q237
Purification:	Affinity-chromatography

Applications

Verified Activity:	Western blot analysis of TRAF6BP expression in HepG2 cell lysate.
Application:	ICC/IF,WB
Recommended	ICC/IF: 1:50-1:200; WB: 1:1000-1:2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Tax1 binding protein 1
Synonyms:	D6Ert404e;CALCOCO3;TXBP151;D6Ert4772e;tax1bp1b;TAX1BP1;T6BP;PRO0105

Research Background

Ubiquitin-binding adapter that participates in inflammatory, antiviral and innate immune processes as well as selective autophagy regulation. Plays a key role in the negative regulation of NF-kappa-B and IRF3 signalings by acting as an adapter for the ubiquitin-editing enzyme A20/TNFAIP3 to bind and inactivate its substrates.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481