

## Anti-Tax1 binding protein 3 Antibody (5F386)

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Clone:	5F386
Purification:	Affinity-chromatography

### Applications

Verified Activity:	Western blot analysis of TAX1BP3 expression in HT-1080 cell lysate.
Application:	FCM,IHC,WB
Recommended	FCM: 1:20-1:100; IHC: 1:100-1:200; WB: 1:1000-1:2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: human Tax1 binding protein 3
Synonyms:	Tax1bp3;TIP1;Tax interaction protein 1;Tax1 binding protein 3;Glutaminase interacting protein 3

### Research Background

May regulate a number of protein-protein interactions by competing for PDZ domain binding sites. Binds CTNNB1 and may thereby act as an inhibitor of the Wnt signaling pathway. Competes with LIN7A for KCNJ4 binding, and thereby promotes KCNJ4 internalization. May play a role in the Rho signaling pathway. May play a role in activation of CDC42 by the viral protein HPV16 E6.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481